

Exercise - 7

October 26, 2023

Objective

Create an iWatch app that reads heart rate data and displays it real-time.

Motivation

You will learn how to use: HealthKit, watch app, .

Must Follow

You must follow the rules below. Otherwise, you will get **50% of your actual score.**

1. Start your XCode project “**Exercise7_LastName_FirstName**” (replace **LastName** with your last name and **FirstName** with your first name).
 - a. **DON'T DO** the following:
 - i. Start with any other project name and change the zip file name later.
 - ii. This will not be accepted at all.
2. You must have to do the exercise compatible with **XCode version 14.1**.

Tips

- **Read the question carefully, then start coding!**
- **Build, Build, and Build**
 - If you add anything on storyboard -> **Build**
 - If you make a reference from storyboard -> **Build**
 - Do not wait until finishing all parts and build.
 - It is easier to debug after each single feature added.

Details

Create a single view iWatch application using Swift as a programming language. Start your XCode project “**Exercise7_LastName_FirstName**” (replace **LastName** with your last name and **FirstName** with your first name).

[5 pts] Design your interface to look like the screenshots [Figures 1-3]

- Design your UI for the Apple Watch
- Single button with the heart emoji, system font, size 50
- Heart Rate value, system font, size 70, regular weight
- “BPM” text, headline font, size 28, bold weight

[2 pts] When user launches the app for the first time it creates permission requests [Figures 4-6].

[1 pt.] After permissions were granted, app is ready to collect HR data (works only for real watch). Simulator will show “0” [Figure 1]

[2 pts] When user taps heart button, emulated outdoor running heart rate data is generated. You might need to restart the app after granting permissions to receive emulated data [Figures 2-3]

- The application should not crash permissions were not granted. Add an appropriate alert for the user in case of any error.

- Make sure to add the appropriate properties to the target properties list:

Privacy – Health Share Usage Description

Privacy – Health Update Usage Description

Privacy – Health Records Usage Description

Bonus

[1 pt.] Create a paired phone app (with a single text label) and transfer heart-rate data to the phone in real time, at least once every 2 sec. (no particular design requirements here, you may add heart to make it similar looking)

[1 pt.] Add data simulation control from the watch (start/stop)

[1 pt.] Add data simulation control from the phone (start/stop)

Submission

Zip XCode project and submit to the blackboard. The name of your zip file will be automatically “**Exercise7_LastName_FirstName.zip**” (**Last Name** is your last name and **FirstName** is your first name). One submission per person.




Figure 1



Figure 2




Figure 3




Figure 4




Figure 5




Figure 6

HINTS

As HealthKit will not be fully covered in the class, here you may find some useful code-snippets. Copy-paste protection activated 😊

1. Perform some function on view appearance:

```
.onAppear(perform: start)
```

2. Request for the authorization:

```
func authorizeHealthKit() {
    let healthKitTypes: Set = [
        HKObjectType.quantityType(forIdentifier:
            HKQuantityTypeIdentifier.heartRate)!]

    healthStore.requestAuthorization(toShare: healthKitTypes, read:
        healthKitTypes) { _, _ in }
}
```

3. Starting query to collect data.

```
private func startHeartRateQuery(quantityTypeIdentifier: HKQuantityTypeIdentifier) {
    // 1 collect data from the device
    let devicePredicate = HKQuery.predicateForObjects(from: [HKDevice.local()])
    // 2 updating handler
    let updateHandler: (HKAnchoredObjectQuery?, [HKSample]?, [HKDeletedObject]?, HKQueryAnchor?, Error?) -> Void = {
        query, samples, deletedObjects, queryAnchor, error in

        guard let samples = samples as? [HKQuantitySample] else {
            return
        }

        self.process(samples, type: quantityTypeIdentifier)

    }

    let query = HKAnchoredObjectQuery(type: HKObjectType.quantityType(forIdentifier: quantityTypeIdentifier)!, predicate:
        devicePredicate, anchor: nil, limit: HKObjectQueryNoLimit, resultsHandler: updateHandler)

    query.updateHandler = updateHandler

    healthStore.execute(query)
}
```

4. Creating workout emulation. That is for the 'heart' button.

```
private func createWorkoutEmulation(){
    // create configuration
    let configuration = HKWorkoutConfiguration()
    configuration.activityType = .running
    configuration.locationType = .outdoor
    // init store, create session and builder variables
    let healthStore = HKHealthStore()
    var session: HKWorkoutSession!
    var builder: HKWorkoutBuilder!
    // create session
    do {
        session = try HKWorkoutSession(healthStore: healthStore, configuration: configuration)
        builder = session.associatedWorkoutBuilder()
    } catch {
        return
    }
    // start session and builder
    session.startActivity(with: Date())
    builder.beginCollection(withStart: Date()) {(success, error) in
    }
}
```

5. Process last heart rate value

```
private func process(_ samples: [HKQuantitySample], type: HKQuantityTypeIdentifier) {
    var lastHeartRate = 0.0

    for sample in samples {
        if type == .heartRate {
            lastHeartRate = sample.quantity.doubleValue(for: heartRateQuantity)
        }
        self.value = Int(lastHeartRate)
    }
}
```

6. Function to run after the view will appear:

```
func start() {
    // authorize monitoring
    authorizeHealthKit()

    // begin collecting heart data
    startHeartRateQuery(quantityTypeIdentifier: .heartRate)
}
```

Again, don't forget that data simulation might require app restart.